


The occlusal force



Dental prescale (Fiji Film, Japan)



The microcapsules are collapsed and the color former leaked from it.



The sheet is bitten by a person.



All occlusal contact

Occluzer FDP-707 (Fuji Film, Japan)

Occlusal force

The occlusal force was calculated after scanning the sheet with an image scanner.



Dental Prescale 50H R type Fuji Film Co., Tokyo, Japan

Bilateral maximal bite force was measured with pressure sensitive sheets.



Maximal Clenching



咬合力の低下

